

Medical Equipment Inspection Service

Anesthesia Machine Service

Scope of Service

- ✓ Inspections performed comply with manufacturer testing specifications
- ✓ Includes periodic replacement of rubber goods, sensors etc.
- ✓ We service Drager, Ohmeda/GE, Dupaco, American Anesthesia, Matrix, Porter, Datascope, and Penlon
- ✓ Servicing includes tasks as noted above in Preventive Maintenance Inspection.
- ✓ Testing to including High and Low Systems pressure and operational testing.
- ✓ Absorber System pressure testing and functionality testing of APL, Inspiratory/ Expiratory valves.
- ✓ Operational testing of Ventilator delivery and monitoring systems.
- ✓ Vaporizer efficacy testing (calibration verification of delivery within manufacturer tolerances).
- ✓ Scavenger System (gas evacuation) inspection, cleaning, and adjustment on anesthesia cart.
- ✓ Complete operational testing of all integrated physiologic monitors and alarm modules.
- ✓ Complete Preventive Maintenance Service and Electrical Safety Inspection

