

Medical Equipment Inspection Service

ELECTRICAL SAFETY/ OPERATIONAL INSPECTION (ESI/OPS) 5 Point Service Scope of Service

California Code of Regulations Title 22, Division. 5, Chapter 1, Article 8, Sec. 70837/70853

Labor for Electrical safety inspections and functional testing will include the following:

- ✓ External physical inspection of the equipment and all of its components.
- ✓ Using an approved electrical safety analyzer to test electrical integrity, power cord continuity, resistance and hazard current leakage of equipment in-patient care areas
- ✓ Perform functional testing, calibration verification and electrical safety testing (ESI) of the unit in accordance with NFPA, TJC, AAMI, OSHA, FDA, and industry standards.
- ✓ Written documentation of service and appropriate equipment labeling.
- ✓ Service summary detailing work performed and applicable Device Safety Alerts

