

Medical Equipment Inspection Service

PERIODIC MAINTENANCE INSPECTION SERVICE (PMs)

Scope of Service

Labor for the periodic service will include the following:

- ✓ Internal and external physical inspection of the equipment and all of its components.
- ✓ Exterior and interior cleaning of debris, filters, hoses and removal of contact oxidation.
- ✓ Lubrication of bearings, motors, pumps, pins, and similar components as necessary.
- ✓ Minor repairs or adjustments requiring up to thirty minutes in duration (labor).
- ✓ Complete functional testing including calibration verification.
- ✓ Electrical Safety Inspection ((ESI/OPS) as noted in ESI 3 point service
- ✓ Appropriate equipment labeling, including measurable output labeling.
- ✓ Complete written documentation of service performed.
- ✓ Service summary detailing work performed and applicable Device Safety Alerts

